

)

Dental plus





Table of Contents...

Restorative : 1-Etch plus

Endodontic : 1-Edta plus solution 2-Edta cream plus 3-Root paste 4-Root pex

Pain management : 1-Candy plus



The purpose of this catalog

This catalog is designed to showcase our full product line. with detailed specifications, packaging options, and value-added features. It also includes exclusive offers and partnership benefits for distributors and clinics.





Catalog Introduction...



- Welcome to Dental Plus, where we redefine the quality and innovation of locally manufactured dental materials with international standards.
- Since our inception, our mission has been clear
- To deliver safe, effective, and scientifically advanced dental products that enhance the experience of both dentists and patients.
- Our vision is to become a leading brand in dental manufacturing across the Middle East and to proudly represent Egyptian-made excellence on a global scale. What sets Dental Plus apart is our commitment to:
- Continuous research and development to ensure product innovation and stability
- Strict in-process and final quality control systems (IPC & QC) at every stage of production
- Reliable after-sales support and technical assistance for professionals and distributors
- Unique product features



Name of product: Etch plus

CSATH HOLS

ETCH PLUS

ETCH PLUS 7

ETCHPUST

TCH PLUS

Code : 12112022

Composition: 37.5% phosphoric acid

Properities :

The etching gel for a strong bond Etch plus is an etching gel containing 37% phosphoric acid. Etch plus is used for enamel etching and dentin conditioning in conjunction with sealants, composite restorations and the adhesive luting of indirect all-ceramic or composite restorations.



Name of product : Etch plus

Advantages

•Balanced handling properties

•Quick, effective application due to consistent material thickness

•Colour can be easily differentiated from the tooth

Special displacer on the plunger ensures economical dispensing of material

•Practical economic price

Area of application

Etch plus is used in the enamel etch or total etch technique for the following indications:

•Composite restorations

Adhesive luting of inlays, crowns, bridges, veneers, adhesive bridges and root canal posts
 Splints

•Luting of orthodontic appliances (e.g. brackets) and bonding of tooth jewellery (e.g. Skyce) •Fissure sealings

DELIVERY FORMS Syringe 3 ml Syringe 5 ml



MG00001

ETCHPLUS

- (LOTI





Name of product: EDTA PLUS SOLUTION

Code : 02022023

Composition: 17% (ethylenediaminetetraacetic acid)

Properities :

EDTA 17% solution with a powerful chelating action for enhanced cleaning and lubrication of the root canals.

Access to tight canal apices Optimal 8,5 pH for optimal performances Reduces friction of endodontic instruments



Name of product: EDTA PLUS SOLUTION

Indications

oEDTA solution 17% (ethylenediaminetetraacetic acid) fi a powerful chelating solution to remove inorganic debris (smear layer) and lubricate the canals in endodontics. oLoosen up calcifications oRemove smear layers oAllow instrumentation (canal shaping) and debridement in a tight/calcified root canal towards the apex.

•Features & benefits oEDTA 17% concentration effectively removes inorganic material and provides possibility to obturate greater number of lateral canals oOptimized chelation provides optimal performance for proper root canal cleaning

DELIVERY FORM : 120 ML BOTTEL





Name of product: Edta plus cream

Code : 04042023

Composition: 17% (ethylenediaminetetraacetic acid) 10% carbameid

Properities :

EDTA cream for simple and effective preparation of the root canals and removal of smear layer and calcification.

Removes inorganic substances

Reaches narrow canals

Reduces friction of endodontic instruments





Name of product: Edta plus cream

Indications

oRoot canal preparation with endodontic instruments oPresence of calcification in the canal oNarrow canals in posterior teeth oSmear layer removal •Features & benefits oReduces friction of endodontic instruments thanks to a lubrication effect for an easy root canal preparation and filing oRemoval of inorganic substances through chelating action for an effective root canal preparation and cleaning

Delivery form : 7 gm syringe





Name of product: Root paste

ROOT PAST

Code : 08082024

Composition: calcium hydroxide , barium sulfate , water

DESCRIPTION:

High pH value for temporary root canal dressing Root paste is a radiopaque, creamy and ready-to-use calcium hydroxide paste for temporary root canal dressing. When temporisation of the root canal between two appointments is required, calcium hydroxide paste is the standard material to use. The product is intended to eliminate remaining microorganisms and prevent reinfection between. The main properties of the calcium hydroxide paste are thanks to its alkaline pH. The usual duration of temporisation treatment is between 1 and 3 weeks. The product features a high pH (above 12). With a radiopacity > 300% of the dentine radiopacity, Opacal is one of the most radiopaque calcium hydroxide paste, ensuring easy identification and efficient follow-up. Root paste is an aqueous non-setting material. This characterisitic induces rapid ionic release and performant action in short treatment time procedurefi Calcium hydroxide pastes: classification and clinical indications). In addition, the paste can be easily removed



Name of product: Root paste

Key benefits

High pH value High radiopacity Creamy consistency Easy and comfortable paste delivery Easy to remove (non-setting material) Rapid ionic release (aqueous composition

Delivery form 3ml syringe (1.5 gm)





Name of product: Root pex

Code : 05052023

Composition: calcium hydroxide , iodoform , silicon oil

Properities Pre-Mixed Calcium Hydroxide paste with Iodoform

INDICATIONS

Temporary root canal treatment Vital pulpotomies in deciduous teeth Apexification and Apexogenesis Stimulation of hard tissue formation Root perforations Root resorption treatment including primary teeth





Name of product: Root pex

FEATURES & BENEFITS

1_Saves time! Pre-mixed and pre-filled
2_Excellent antibacterial and bacteriostatic properties
3_Highly radiopaque
4_Fast and easy application
5_Pharmaceutical grade iodoform
6_Complete coating of canal walls
7_Does not harden, so Ca(OH)2 remains active

Delivery form : syringe 3ml (1.5gm)







Name of product: candy plus

Composition: caine , geling agent

Indications

o Topical analgesia of oral mucosa before anaesthetic injection
• Temporary Relief of pain due to minor irritation or injury of the oral mucosa and gums.
• Temporary relief of pain associated with ,denture or orthodontic appliances
• Gingival Application before Scaling.

•Features & benefits

o Fast and effective Topical anaesthetic gel o Easy to use gel formulation with Dye-Free flavouring and no bitter taste.

Delivery form : 30gms jar





Dental plus

01006399095
 01091259027

DENTAL

DENTALPLUSZI@SGMAIL.COM

HTTPS://www.instagram.com/dentalsplus/

0

HTTPS://WWW.TIKTOK.COM//IDENTALPEUSEG

ð



🧾 الخاتمة والشكر

نشكركم على وقتكم وامتمامكم بكتالوج منتجات Dental plus

ونأمل أن نكون قد قدمنا لكم صورة واضحة عن جودة منتجاتنا وتميزنا في مجال مستلزمات طب الأسنان

فريقنا دانمًا على استعداد للإجابة عن استفساراتكم، وتقديم الدعم الفني، وتوفير أي متطلبات خاصة أو عروض تعاونية



🧾 لا تترددوا في التواصل معنا – سواءً لطلب منتج، أو استشارة، أو زيارة مصنعنا يسعدنا دانقا استقبالكم والتعاون معكم معًا نصنع مستقبلًا أكثر تميزًا في طب الأسنان نتشرف بخدمتكم ونضمن لكم جودة المنتجات وأفضل الأسعار





Thank You & Get in Touch

We truly appreciate your time exploring our catalog and considering Dental Plus as your trusted partner in dental care solutions.

Our team is always ready to assist you with product details, custom requests, or partnership opportunities.

Feel free to reach out fl whether you're looking to place an order, ask a question, or visit our facility, we,,d be honored to welcome you.

Together, let, s build the future of dental excellence.

